


## Guide – Listen to music on your hearing device from your phone

### Before you start to pair your hearing device

- You may have to download an app from the company that makes your hearing aid.
- Your audiologist will tell you which app and how to download.
- Also, make sure that your Bluetooth is turned on in Settings.

### Enable hearing device pairing

- Go to settings.
- Scroll down to Accessibility.
- Tap accessibility.
- Scroll down to until you see Hearing Devices.
- Tap.
- You can see if the iPhone is paired with the hearing aid at the top of the menu.
- If it is not paired (or loses the Bluetooth connection), you may have to take a few steps. First, take off your hearing aid and turn it off for three seconds (Starkey devices).
- Then, open the hearing device app and follow instructions to pair your hearing aid again.



## Change color contrast

- Once paired, go back into your Hearing Device settings on the iPhone.
- Tap on Audio Routing.
- You can then choose your settings for call audio (phone calls) or media (music, videos, and podcasts).
- If you don't want these to go directly to your hearing aid, select "Never".
- If you do want to hear directly into your hearing aids, select "Automatic" or "Always".

